

### **Before the first game:**

1. Sit on the chair and turn on the computer
2. Put on the badge with arduino
3. On the right hand:
  - 3.1. Put on GSR on the right hand (ring and little finger)
  - 3.2. Put on IMU
  - 3.3. Put on EMG electrodes
4. On the left hand:
  - 4.1. Put on IMU
  - 4.2. Put on GSR electrodes
5. Log in to the game
6. Synchronise the time on the computer
7. Set the permanent monitor position
8. Check if webcam can record face expressions
9. Calibrate eye tracker

### **Before each game:**

1. Start the input logger
2. Start screen recording
3. Start webcam recording
4. Start eye tracker recordings
5. Start Sensors recordings
6. Search for a game
7. Write down the exact time when the game starts

### **After each game:**

1. Save replays
2. Turn off input logger
3. Turn off screen recording
4. Turn off webcam recording
5. Turn off eye tracker recordings
6. Turn off sensors recordings

### **After the last game:**

1. Ask each player to fill in the form for Amazon voucher
- 2.